

# **World Leading Supercontinuum to 4.75 μm**

Peter M. Moselund, NKT Photonics

[pmm@nktphotonics.com](mailto:pmm@nktphotonics.com)



Crystal Fibre • *aeroPULSE* • Koheras • SuperK • LIOS • Fianium

**NKT** Photonics

# What is supercontinuum

- Broadband light source – like a lamp or globar
  - High brightness – like a laser
  - Spatially coherent – like a laser
  - Single mode fiber output – like a laser
  - Pulsed light source – like a laser
- **This is why they are nicknamed “white lasers”**

Mid-IR Supercontinuum

- Constructed from telecom components
- Covers the entire spectrum



# Supercontinuum – How?



- |                                                                                                                                                                       |                                                                                                                                                                                                   |                                                                                                                                                                                   |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"><li>➤ ps-ns</li><li>➤ single pulse -100MHz</li><li>➤ kW peak power</li><li>➤ Single mode</li><li>➤ Preferably fiber laser</li></ul> | <ul style="list-style-type: none"><li>➤ 3-10 μm dia. core</li><li>➤ Length: some meters</li><li>➤ Silica, fluoride or chalcogenide</li><li>➤ Single mode</li><li>➤ Anomalous dispersion</li></ul> | <ul style="list-style-type: none"><li>➤ Ultra broadband</li><li>➤ All colors in one beam</li><li>➤ Single mode</li><li>➤ Short pulsed</li><li>➤ Collimated fiber output</li></ul> |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

# Why Supercontinuum


## Positive

- High brightness
- Nothing except synchrotrons can give wider bandwidth
- Great for access to multiple absorption lines
- Very flexible system
- Inherently fiber delivered
- Always pulsed




## Negative

- Spreading the spectrum means less light in each line
- Filtering has to be done on detection side
- Higher power consumption



# Supercontinuum's colors

Established Engineering Department, Aalto University, Finland




## Setting the limit for ZBLAN

ZBLAN is the most mature soft glass

The limit for ZBLAN SC is often named  
at  $4.3\text{ }\mu\text{m}$

- In MINERVA NKT has shown that ZBLAN can deliver SC up to  $4.75\text{ }\mu\text{m}$  even in a compact source

1200 mW  
350 kHz  
3ns  
 $1.7\text{-}4.75\text{ }\mu\text{m}$



**NKT** Photonics  
the power of light

# Mid-IR Source generations



Crystal Fibre • *aeroPULSE* • Koheras • SuperK • LIOS • Fianium




# SuperK MIR



- 2.5 MHz
- 1800 – 4200 nm
- Total power > 450 mW
- Single mode fiber delivery
- $M^2 < 1.1$
- Fiber delivered collimated output
- Room temperature operation
- DEMO source is available now


# Power Density curve



Crystal Fibre • *aeroPULSE* • Koheras • SuperK • LIOS • Fianium

**NKT** Photonics

# Spectral tailoring



Crystal Fibre • *aeroPULSE* • Koheras • SuperK • LIOS • Fianium

**NKT** Photonics


# Scaling up the Power

11

- Total power: 4.8 W
- 60% of power >2.5 μm
- Power >3.5μm: 1.34 W
- Record for power >3.5 μm, at the time



Crystal Fibre • aeroLASE • Koheras • SuperK



**NKT** Photonics  
the power of light

# Filtering


Mid-IR SC has been tested with a number of filtering options

- Monochromators
- Fourier transform spectrometers (FTIRs)
- Acousto optical filters (AOTFs)
- Fabry-Perot micro spectrometers (FPMS)
- Up-conversion detectors



# Reliability

- Working with soft glass fiber is challenging.
- In the first 2-3 years of our development this gave reliability problems.
- Now we have demonstrated over 2 years 24-7 operation
- Over 13,000 hours runtime completely service free
- Operation tested at 0-40 °C




# Applications: Storage atmosphere optimization

- SC driven multi-gas spectroscopy to control storage and ripening of potatoes and fruit.
- QCap, InterReg NWE project, will run until December 2019



- Aims to detect gasses such as:  
CO<sub>2</sub>, O<sub>2</sub>, Ethylene, Ethanol, Ethane, Methanol,  
Acetone, Acetaldehyde, Ethylacetate



# Applications: Drone based gas measurement

- SC driven multigas spectroscopy on a drone for pollution monitoring
- H2020 project, FLAIR will run until December 2019
- Target: 8L, 7kg total system
- NKT will develop world's smallest Supercontinuum source



**FLAIR**  
www.h2020flair.eu



Crystal Fibre • *aeroPULSE* • Koheras • SuperK • LIOS • Fianium

**NKT** Photonics

## Applications in general

- Many diverse scientific applications
- And a few key areas
  - Defence
  - mid-IR tissue/cell analysis
  - Non destructive testing
  - Multi-Gas Spectroscopy

## Conclusion

- Supercontinuum sources can provide high brightness over a wider spectrum than any other portable source
- Mid-IR supercontinuum is now becoming a well proven technology
- Sources are commercially available also in mid-IR
- Applications are expected to grow rapidly in the coming years



# Questions?

Part of this work has been supported by the European Commission under the Seventh Framework Programme (FP7) under grant agreement No 317803  
Other parts of this work has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement No 732968

Crystal Fibre • *aeroPULSE* • Koheras • SuperK • LIOS • Fianium

**NKT** Photonics